

Product Datasheet

Human Carnitine-acylcarnitine translocase (CACT) ELISA Kit (orb1656189)

Description	Human Carnitine-acylcarnitine translocase (CACT) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625-40 ng/mL. Sensitivity: 0.235 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.