

Product Datasheet

Human Secreted Ly-6/uPAR-related protein 1 (SLURP1) ELISA Kit (orb1656246)

Description	Human Secreted Ly-6/uPAR-related protein 1 (SLURP1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.058 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.