

Product Datasheet

Fish Superoxide Dismutase (SOD) ELISA Kit (orb1656259)

Description	Fish Superoxide Dismutase (SOD) ELISA Kit
Reactivity	Fish
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 4-140 IU/mL. Sensitivity: 1.8 IU/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.