

Product Datasheet

GALNT14 antibody (orb165635)

Description	Rabbit polyclonal antibody to GALNT14 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GALNT14
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	64321
Uniprot ID	Q96FL9
NCBI	NP_001240756.1 , NP_001240755.1 , NP_078848.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

