

Product Datasheet

KBTBD7 antibody (orb165638)

Description	Rabbit polyclonal antibody to KBTBD7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KBTBD7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	77163
Uniprot ID	Q8WVZ9
NCBI	NP_115514.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

