

Product Datasheet

Human Homeobox protein TGIF2LY (TGIF2LY) ELISA Kit (orb1656381)

Description	Human Homeobox protein TGIF2LY (TGIF2LY) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.