

Product Datasheet

Mouse Transforming growth factor β 2 (TGF- β 2) ELISA Kit (orb1656412)

Description Mouse Transforming growth factor
 β 2 (TGF- β 2) ELISA Kit

Reactivity Mouse

Tested ELISA

Applications

Note For research use only

Application notes Detect Range: Request
Information. Sensitivity: Request
Information. Assay Time: 1-4.5 h.
Detection wavelength: 450 nm.
Sample Volume: 1-200 μ L. Sample
Type: Serum, Plasma, Other
biological fluids.

Expiration Date Please enquire.