

Product Datasheet

CCNJ antibody (orb165654)

Description	Rabbit polyclonal antibody to CCNJ
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CCNJ
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42581
Uniprot ID	Q5T5M9
NCBI	NP_001127847.1 , NP_001127848.1 , NP_061957.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

