

Product Datasheet

DNAJB5 antibody (orb165656)

Description	Rabbit polyclonal antibody to DNAJB5 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	DNAJB5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	39133
Uniprot ID	O75953
NCBI	NP_001128476.2 , NP_036398.3 , NP_001128477.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

