

Product Datasheet

KIRREL2 antibody (orb165664)

Description	Rabbit polyclonal antibody to KIRREL2 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KIRREL2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	75093
Uniprot ID	Q6UWL6
NCBI	NP_954649.2 , NP_115499.4 , NP_954648.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

