

Product Datasheet

DNAJC10 antibody (orb165671)

Description	Rabbit polyclonal antibody to DNAJC10 (N-term).
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	DNAJC10
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	91080
Uniprot ID	Q8IXB1
NCBI	NP_061854.1 , NP_001258510.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

