

Product Datasheet

USP39 antibody (orb165673)

Description	Rabbit polyclonal antibody to USP39
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	USP39
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	65381
Uniprot ID	Q53GS9
NCBI	NP_001243656.1 , NP_001243654.1 , NP_001243657.1 , NP_006581.2 , NP_001243655.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

