

Product Datasheet

SSBP2 antibody (orb165674)

Description	Rabbit polyclonal antibody to SSBP2 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SSBP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37828
Uniprot ID	P81877
NCBI	NP_001243662.1, NP_036578.2, NP_001243661.1, NP_001243663.1, NP_001243665.1, NP_001243664.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

