

Product Datasheet

Human Signal recognition particle 54 kDa protein (SRP54) ELISA Kit (orb1656819)

Description Human Signal recognition particle 54 kDa protein (SRP54) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.