

Product Datasheet

Parva antibody (orb165685)

Description	Rabbit polyclonal antibody to Parva (N-term).
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Parva
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42330
Uniprot ID	Q9EPC1
NCBI	NP_065631.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

