

Product Datasheet

HAS3 antibody (orb165688)

Description	Rabbit polyclonal antibody to HAS3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	HAS3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62998
Uniprot ID	O00219
NCBI	NP_005320.2 , NP_619515.1 , NP_001186209.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

