

Product Datasheet

FOXP3 antibody (orb165697)

Description Rabbit polyclonal antibody to FOXP3 (C-term S418).

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target FOXP3

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 47244

Uniprot ID [Q9BZS1](#)

NCBI [NP_054728.2](#), [NP_001107849.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

