

## Product Datasheet

# Sheep Adrenocorticotrophic hormone (ACTH) ELISA Kit (orb1657197)

|                            |                                                                                                                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Sheep Adrenocorticotrophic hormone (ACTH) ELISA Kit                                                                                                                                        |
| <b>Reactivity</b>          | Sheep                                                                                                                                                                                      |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                      |
| <b>Note</b>                | For research use only                                                                                                                                                                      |
| <b>Application notes</b>   | Detect Range: 1.56-100 pg/mL.<br>Sensitivity: 0.72 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                            |