

Product Datasheet

Mouse Tenascin (TNC) ELISA Kit (orb1657200)

Description	Mouse Tenascin (TNC) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.78-50 ng/mL. Sensitivity: 0.195 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.