

Product Datasheet

Cat Tumor necrosis factor (TNF/TNFA/TNFSF2) ELISA Kit (orb1657204)

Description	Cat Tumor necrosis factor (TNF/TNFA/TNFSF2) ELISA Kit
Reactivity	Feline
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.8-500 pg/mL. Sensitivity: 3.1 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.