

Product Datasheet

Human Tumor necrosis factor receptor superfamily member 21 (TNFRSF21) ELISA Kit (orb1657255)

Description Human Tumor necrosis factor receptor superfamily member 21 (TNFRSF21) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL. Sensitivity: 0.056 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.