

Product Datasheet

TPD52L1 antibody (orb165727)

Description	Rabbit polyclonal antibody to TPD52L1 (N-term).
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TPD52L1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22449
Uniprot ID	Q16890
NCBI	NP_001003397.1 , NP_001003395.1 , NP_001003396.1 , NP_003278.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

