

Product Datasheet

Fish Tumor Necrosis Factor α (TNF- α) ELISA Kit (orb1657281)

Description	Fish Tumor Necrosis Factor α (TNF- α) ELISA Kit
Reactivity	Fish
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.8-500 pg/mL. Sensitivity: 3.7 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.