

Product Datasheet

Cat Troponin I, cardiac muscle (TNNI3) ELISA Kit (orb1657319)

Description	Cat Troponin I, cardiac muscle (TNNI3) ELISA Kit
Reactivity	Feline
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 2.34-150 ng/mL. Sensitivity: 0.65 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.