



Product Datasheet

Rat Stra6 antibody (orb165732)

Description	Rabbit polyclonal antibody to Rat Stra6 (C-term).
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Stra6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	73763
Uniprot ID	Q4QR83
NCBI	NP_001025095.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat heart tissue...