

Product Datasheet

Plant Vitamin A (VA) ELISA Kit (orb1657561)

Description	Plant Vitamin A (VA) ELISA Kit
Reactivity	Plant
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625-10 ng/mL. Sensitivity: 0.467 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.