

Product Datasheet

Rabbit Von Willebrand Factor (VWF) ELISA Kit (orb1657627)

Description	Rabbit Von Willebrand Factor (VWF) ELISA Kit
Reactivity	Rabbit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 12.35-1000 ng/mL. Sensitivity: 5.25 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.