

Product Datasheet

Rat BMP-9 ELISA Kit (orb1657699)

Description	Rat BMP-9 ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: Gdf2. Gene_full_name: growth differentiation factor 2. Protein_name: Growth differentiation factor 2. Sample Type: bone tissue cell culture supernates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat.
Expiration Date	Please enquire.