

Product Datasheet

PPP1R9A antibody (orb165773)

Description	Rabbit polyclonal antibody to PPP1R9A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PPP1R9A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	123342
Uniprot ID	Q9ULJ8
NCBI	NP_001159632.1 , NP_001159634.1 , NP_001159635.1 , NP_060120.2 , NP_001159633.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

