

Product Datasheet

Rat TNFSF11/RANKL ELISA Kit (orb1657731)

Description	Rat TNFSF11/RANKL ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6 pg/mL-1000 pg/mL. Gene_name: Tnfsf11. Gene_full_name: Tumor necrosis factor ligand superfamily member 11 . Protein_name: Tumor necrosis factor ligand superfamily member 11. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Tissue_specificity: Highly expressed in thymus and bone tissues.
Expiration Date	Please enquire.