

Product Datasheet

IL15RA antibody (orb165782)

Description	Rabbit polyclonal antibody to IL15RA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	IL15RA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	28233
Uniprot ID	Q13261
NCBI	NP_001230468.1 , NP_002180.1 , NP_001243694.1 , NP_751950.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

