

Product Datasheet

Mouse TNFSF18/GITRL ELISA Kit (orb1657877)

Description	Mouse TNFSF18/GITRL ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: TNFSF18. Gene_full_name: Tumor necrosis factor ligand superfamily member 18. Protein_name: Tumor necrosis factor ligand superfamily member 18. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity cytokines necrosis factors tnf superfamily tags & cell markers cell type markers epi endo-thelial signal transduction growth factors hormones tnf cancer growth factors cardiovascular vasculature endothelium. Tissue_specificity: Detected in immature and mature dendriticcells and in macrophages (at protein level). Detected in spleen lung heart thymus monocytes macrophages B-cells and dendriticcells.
Expiration Date	Please enquire.