

Product Datasheet

Mouse CD22 ELISA Kit (orb1657882)

Description	Mouse CD22 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 156 pg/mL-10000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: CD22. Gene_full_name: B-cell receptor CD22. Protein_name: B-cell receptor CD22. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology adaptive immunity b cells cd cancer tumor immunology tumor-associated antigens tumor biomarkers tumor antigens stem cells hematopoietic progenitors lymphoid b lymphocytic lineage. Sequence_similarities: Belongs to the immunoglobulin superfamily. SIGLEC (sialic acid binding Ig-like lectin) family. Tissue_specificity: B-lymphocytes.
Expiration Date	Please enquire.