

Product Datasheet

Rat ADAM12 ELISA Kit (orb1657890)

Description	Rat ADAM12 ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 156 pg/mL-10000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: ADAM12. Gene_full_name: ADAM metallopeptidase domain 12. Protein_name: Disintegrin and metalloproteinase domain-containing protein 12. Sample Type: cell culture supernates serum and plasma (heparin). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: signal transduction cytoskeleton ecm cell adhesion integrins alpha beta extracellular matrix ecm enzymes adam protein family cancer invasion microenvironment cell biology proteolysis ubiquitin proteolytic enzymes metalloprotease adams. Sequence_similarities: Contains 1 disintegrin domain.
Expiration Date	Please enquire.