

Product Datasheet

Human CRLF2/TSLP R ELISA Kit (orb1657893)

Description	Human CRLF2/TSLP R ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 156 pg/mL-10000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: CRLF2. Gene_full_name: Cytokine receptor-like factor 2. Protein_name: Cytokine receptor-like factor 2. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity cytokines. Sequence_similarities: Belongs to the type I cytokine receptor family. Type 5 subfamily. Tissue_specificity: Expressed in heart skeletal muscle kidneyand adult and fetal liver. Primarily expressed in dendrites andmonocytes. Weakly expressed in T-cells.
Expiration Date	Please enquire.