

Product Datasheet

FRS2 antibody (orb165791)

Description	Rabbit polyclonal antibody to FRS2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	FRS2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57029
Uniprot ID	Q8WU20
NCBI	NP_001036020.1 , NP_006645.3 , NP_001265286.1 , NP_001265283.1 , NP_001265285.1 , NP_001265284.1 , NP_001265280.1 , NP_001265282.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

