

Product Datasheet

Mouse JAM-A/F11R ELISA Kit (orb1657925)

Description	Mouse JAM-A/F11R ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: F11R. Gene_full_name: Junctional adhesion molecule A. Protein_name: Junctional adhesion molecule A. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: cardiovascular blood platelets signal transduction cytoskeleton ecm cell adhesion tight junctions atherosclerosis vascular inflammation leukocyte recruitment cell junction molecules. Sequence_similarities: Belongs to the immunoglobulin superfamily.
Expiration Date	Please enquire.