

Product Datasheet

Human TNFRSF13B/TACI ELISA Kit (orb1658010)

Description	Human TNFRSF13B/TACI ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 93.8 pg/mL-6000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: TNFRSF13B. Gene_full_name: Tumor necrosis factor receptor superfamily member 13B. Protein_name: Tumor necrosis factor receptor superfamily member 13B. Sample Type: cell culture supernates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: Nothing Found. Sequence_similarities: Contains 2 TNFR-Cys repeats. Tissue_specificity: Highly expressed in spleen thymus smallintestine and peripheral blood leukocytes. Expressed in resting B-cells and activated T-cells but not in resting T-cells.
Expiration Date	Please enquire.