

Product Datasheet

Mouse CSF1R/M-CSFR ELISA Kit (orb1658113)

Description	Mouse CSF1R/M-CSFR ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5 pg/mL-4000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: CSF1R. Gene_full_name: Macrophage colony-stimulating factor 1 receptor. Protein_name: Macrophage colony-stimulating factor 1 receptor. Sample Type: cell culture supernates cell lysates tissue homogenates serum or plasma (heparin EDTA). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: cardiovascular angiogenesis cytokines immunology innate immunity macrophage inflamm. csfs signal transduction protein phosphorylation ser thr kinases kinases. Sequence_similarities: Belongs to the protein kinase superfamily. Tyr protein kinase family. CSF-1 PDGF receptor subfamily. Tissue_specificity: Widely expressed.
Expiration Date	Please enquire.