

Product Datasheet

Mouse MBL2/MPB-C ELISA Kit (orb1658115)

Description	Mouse MBL2/MPB-C ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 312 pg/mL-20 000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: MBL2. Gene_full_name: Mannose-binding protein C. Protein_name: Mannose-binding protein C. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: MBL2. Sequence_similarities: Contains 1 C-type lectin domain.
Expiration Date	Please enquire.