

Product Datasheet

Mouse HAI-2/SPINT2 ELISA Kit (orb1658134)

Description	Mouse HAI-2/SPINT2 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2 pg/mL-2000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: Spint2. Gene_full_name: Kunitz-type protease inhibitor 2. Protein_name: Kunitz-type protease inhibitor 2. Sample Type: cell culture supernates serum plasma (heparin EDTA) cell lysates and tissue homogenates. Tmb_incubation: 15-20 min. Research_category: SPINT2. Tissue_specificity: Isoform 2 is more predominantly expressed than isoform 1.
Expiration Date	Please enquire.