

Product Datasheet

Mouse SPINT1/HAI-1 ELISA Kit (orb1658136)

Description	Mouse SPINT1/HAI-1 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5 pg/mL-4000 pg/mL. Sensitivity: 12 pg/mL. Gene_name: Spint1. Gene_full_name: Kunitz-type protease inhibitor 1. Protein_name: Kunitz-type protease inhibitor 1. Sample Type: cell culture supernates serum and plasma (heparin). Tmb_incubation: 15-20 min. Research_category: cell biology proteolysis ubiquitin protease inhibitors protease inhibitors.
Expiration Date	Please enquire.