

Product Datasheet

PPIAL4A antibody (orb165816)

Description	Rabbit polyclonal antibody to PPIAL4A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PPIAL4A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	18182
Uniprot ID	Q9Y536
NCBI	NP_839944.1 , NP_001129261.1 , NP_001137355.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

