

Product Datasheet

Human CXCL4/PF4 ELISA Kit (orb1658166)

Description	Human CXCL4/PF4 ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 156 pg/mL-10000 pg/mL. Sensitivity: 15 pg/mL. Gene_name: PF4. Gene_full_name: Platelet factor 4. Protein_name: Platelet factor 4. Sample Type: cell culture supernates cell lysates tissue homogenates serum and plasma (heparin EDTA). Tmb_incubation: 15-20 min. Research_category: cardiovascular angiogenesis inhibitors immunology innate immunity macrophage inflamm. blood platelets chemokines alpha chemokines (cxc) cancer invasion microenvironment angiogenic inhibitory factors atherosclerosis vascular inflammation leukocyte recruitment kits lysates other kits elisa kits atherosclerotic proteins elisa kits.
Expiration Date	Please enquire.