

Product Datasheet

Mouse IL7R/CD127 ELISA Kit (orb1658173)

Description	Mouse IL7R/CD127 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 312 pg/mL-20 000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: IL7R. Gene_full_name: Interleukin-7 receptor subunit alpha. Protein_name: Interleukin-7 receptor subunit alpha. Sample Type: cell culture supernates and serum. Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity cytokines interleukins stem cells hematopoietic progenitors lymphoid b lymphocytic lineage t lymphocytic lineage myeloid monocytic lineage. Sequence_similarities: Belongs to the type I cytokine receptor family. Type 4 subfamily. Tissue_specificity: Spleen thymus and fetal liver.
Expiration Date	Please enquire.