

Product Datasheet

Rat B7-1/CD80 ELISA Kit (orb1658177)

Description	Rat B7-1/CD80 ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5 pg/mL-4000 pg/mL. Sensitivity: 2 pg/mL. Gene_name: CD80. Gene_full_name: CD80 molecule. Protein_name: T-lymphocyte activation antigen CD80. Sample Type: cell culture supernates and serum. Tmb_incubation: 15-20 min. Capture_antibody: monoclonal antibody from mouse. Detection_antibody: polyclonal antibody from goat. Research_category: immunology adaptive immunity b cells cd stem cells hematopoietic progenitors lymphoid b lymphocytic lineage myeloid dendritic cell lineage monocytic lineage.
Expiration Date	Please enquire.