

Product Datasheet

PRAMEF16 antibody (orb165819)

Description	Rabbit polyclonal antibody to PRAMEF16 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	This PRAMEF16 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 446-472 amino acids from the C-terminal region of human PRAMEF16.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q5VWM1
NCBI	NP_001038945.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

