

Product Datasheet

Mouse CD30/TNFRSF8 ELISA Kit (orb1658225)

Description	Mouse CD30/TNFRSF8 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6 pg/mL-1000 pg/mL. Sensitivity: 10 pg/mL. Gene_name: TNFRSF8. Gene_full_name: Tumor necrosis factor receptor superfamily member 8. Protein_name: Tumor necrosis factor receptor superfamily member 8. Sample Type: cell culture supernates and serum. Tmb_incubation: 25-30 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology cell type markers cd non-lineage innate immunity cytokines tnf superfamily cancer tumor immunology cd markers tumor biomarkers receptors. Tissue_specificity: Detected in thymus and in activatedsplenocytes.
Expiration Date	Please enquire.