

Product Datasheet

Mouse MIP-1Alpha/CCL3 ELISA Kit (orb1658280)

Description	Mouse MIP-1Alpha/CCL3 ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.8 pg/mL-500 pg/mL. Sensitivity: 10 pg/mL. Gene_name: CCL3. Gene_full_name: C-C motif chemokine 3. Protein_name: C-C motif chemokine 3. Sample Type: cell culture supernates serum and plasma (heparin EDTA citrate). Tmb_incubation: 20-25 min. Capture_antibody: monoclonal antibody from rat. Detection_antibody: polyclonal antibody from goat. Research_category: immunology innate immunity macrophage inflamm. chemokines beta chemokines (cc) kits lysates other kits elisa kits blood cell antigens elisa kits immune system diseases antiviral signaling hiv-related. Tissue_specificity: Expressed in lung spleen and pancreas.
Expiration Date	Please enquire.