

Product Datasheet

VSIG4 antibody (orb165831)

Description	Rabbit polyclonal antibody to VSIG4 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	VSIG4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43987
Uniprot ID	Q9Y279
NCBI	NP_001171759.1 , NP_001093901.1 , NP_009199.1 , NP_001244332.1 , NP_001171760.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

